

ACCESSORIES

FOR MEDICAL DIAGNOSTIC EQUIPMENT



ACCESSORIES FOR MEDICAL DIAGNOSTIC EQUIPMENT

Neurosoft Company has been manufacturing and supplying the medical facilities of the world with the complex equipment for neurophysiology and electrodiagnostics for more than 20 years.

We guarantee our customers the high quality of our devices, constantly updated software and after-sale service during the whole warranty period.

However to ensure the optimal performance, the device shall be equipped with the high-quality and reliable accessories and consumables. That is why we manufacture the most part of the accessories for our devices ourselves using the components from Russia, Germany, France, the USA and Italy. Each produced item is carefully checked and tested together with the supplied device.

Besides we purchase and supply the qualitative video and computer equipment compatible with our devices.

You have the Neurosoft catalogue in your hands. It describes in details the supplied items, their specifications and also the order code. You can just select the required item, inform the manager of sales department using e-mail address indicated on the last page of the catalogue and be sure that you get exactly what you need.

* We reserve the right to improve the design and specifications of the described products at any time without previous warning.

CONTENTS

● EEG Accessories

Basic Accessories	5
PSG Accessories	6
Electrode Systems for Long-term EEG Monitoring	9
Video EEG Accessories	11
Stimulators	13
Miscellaneous	14

● EMG Accessories

Basic Accessories	16
Stimulators	18
Miscellaneous	21

● IONM Accessories

Basic Accessories	24
Electrodes	25
Miscellaneous	27

● ABR and OAE Accessories

Basic Accessories	29
Stimulators	30
Miscellaneous	32

● ERG Accessories

Basic Accessories	34
Stimulators	35
Miscellaneous	35

● Accessories for Magnetic Stimulators

Coils for Magnetic Stimulators	37
Basic Accessories	38

● ECG Accessories

Basic Accessories	41
Miscellaneous	43

● Accessories for Spirometers

Basic Accessories	45
-------------------------	----

● Stands, Arms, Holders and Trolleys, Pastes and Gels

Stands, Arms, Holders	47
Trolleys	50
Pastes, Gels	51

EEG ACCESSORIES



BASIC ACCESSORIES

Bridge EEG electrode



ORDER CODE: NS015106.007

Ear EEG electrode



ORDER CODE: NS015106.015

Cable for bridge or ear EEG electrode



ORDER CODE	COLOUR
NS990103.036-01.10	○
NS990103.036-03.10	●
NS990103.036-04.10	●

LENGTH: 1 m

Ear cup EEG electrode with cable



ORDER CODE: NS015106.018

Cup EEG electrode with cable



ORDER CODE	COLOUR
NS990106.072-01.10	○
NS990106.072-03.10	●

LENGTH: 1 m

EEG helmet (adult)



ORDER CODE	SIZE	COLOUR
NS023211.007-00	42-48	●
NS023211.007-01	48-54	●
NS023211.007-02	54-62	●

EEG helmet (pediatric)



ORDER CODE	SIZE
NS023211.002-01	32-38
NS023211.002-02	34-42

Cup electrode with cable



ORDER CODE	COLOUR
NS990106.028-01.10	○
NS990106.028-02.10	●
NS990106.028-03.10	●
NS990106.028-04.10	●
NS990106.028-05.10	●
NS990106.028-06.10	●
NS990106.028-09.10	●
NS990106.028-10.10	●

LENGTH: 1 m

Reusable cup electrode with cable C1234 (silver)



ORDER CODE	DIAMETER	LENGTH
NS990106.049-001	6 mm	1 m
NS990106.049-002	6 mm	1.5 m
NS990106.049-003	6 mm	2 m
NS990106.049-004	10 mm	1 m
NS990106.049-005	10 mm	1.5 m
NS990106.049-006	10 mm	2 m

OTHER CUP COATING:

- silver-chloride
- gold

Disposable subdermal single needle electrode with cable S50716



ORDER CODE: NS990106.043-001

LENGTH: 1 m

Disposable subdermal single needle electrode with cable S50716 (twisted pair)



ORDER CODE: NS990106.043-005

LENGTH: 2 m

Y-adapter



ORDER CODE: NS006103.019

Disposable ECG electrode with cable F9058N (white)



ORDER CODE: NS226934

SIZE: 23×30 mm
QTY: 3 pcs.

Disposable ECG electrode with coloured cable F9058NRM



ORDER CODE: NS103096

SIZE: 22 mm
QTY: 3 pcs.

PSG ACCESSORIES

Neuron-Spectrum-PU5 patient unit



ORDER CODE: NS055998.001

KIT CONTAINS:

- Neuron-Spectrum-PU5 patient unit (NS055201.004)
- Cable for patient unit connection
- Neck strap (NS990211.007)
- Technical manual

FOR THE FOLLOWING DEVICES:

- Neuron-Spectrum-1, 2, 3, 4
- Neuron-Spectrum-4/P
- Neuron-Spectrum-4/EPM
- Neuron-Spectrum-5, 5/S

Cup electrode with cable (for limb movement recording)



ORDER CODE: NS990106.028-01.30

LENGTH: 3 m

Cable for disposable electrode with button connector



ORDER CODE	LENGTH	COLOUR
NS990103.037-03.10	1 m	●
NS990103.037-03.15	1.5 m	●
NS990103.037-03.25	2.5 m	●
NS990103.037-04.10	1 m	●
NS990103.037-04.15	1.5 m	●
NS990103.037-04.25	2.5 m	●

CONNECTOR: touch-proof

Airflow sensor ATS-1-20



ORDER CODE: NS990356.004

CONNECTOR: 3.5 stereo
LENGTH: 2 m

Airflow pressure sensor APS-1



ORDER CODE: NS990998.022

KIT CONTAINS:

- Airflow pressure sensor APS-1
- Nasal cannula
- Adapter for nasal cannula MTL230-9

CONNECTOR: touch-proof
LENGTH: 2 m

Airflow pressure sensor APS-2



ORDER CODE: NS990998.022-01

KIT CONTAINS:

- Airflow pressure sensor APS-2
- Nasal cannula
- Adapter for nasal cannula MTL230-9

CONNECTOR: 3.5 mm stereo
LENGTH: 2 m

Nasal cannula



ORDER CODE: NS104216


LENGTH: 2 m

Adapter for nasal cannula MTL230-9



ORDER CODE: NS222562

Snoring sensor



ORDER CODE	LENGTH	MODEL
NS990356.006	2 m	SS-1
NS990356.006-01	0.8 m	SS-1-08 (for Neuron-Spectrum-AM)

Chest or abdominal respiratory effort sensor



ORDER CODE	LENGTH	MODEL
NS990998.020	2 m	ES-4-22
NS990998.020-02	0.5 m	ES-4-05 (for Neuron-Spectrum-AM)

CONNECTOR: touch-proof

Body position sensor PS-1



ORDER CODE	CONFIGURATION
NS990998.012-01	with belts and adapter
NS990356.003	without belts and adapter

**Adapter
2 touch-proof → 3.5 stereo**




ORDER CODE: NS990103.024

Reusable finger sensor for SpO2 recording (adult)




ORDER CODE: NS208851

Reusable finger sensor for SpO2 recording (pediatric)



ORDER CODE: NS216409

Disposable finger sensor for SpO2 recording



ORDER CODE: NS228039

Reusable limb clamp ECG electrode



ORDER CODE	MODEL	SIZE
NS103212	F9024SSC	adult
NS101823	F9023	pediatric

QTY: 4 pcs.

3-lead ECG cable



ORDER CODE: NS007103.016

External unit for CO2 recording



ORDER CODE: NS051999.001



KIT CONTAINS:

- External unit for CO2 recording (NS051201.001)
- Anaesthetic face mask (NS102257)
- Technical manual

External unit for SpO2 recording



ORDER CODE: NS045201.001

ELECTRODE SYSTEMS FOR LONG-TERM EEG MONITORING

Electrode system Electro-Cap



ORDER CODE	SIZE	COLOUR
NS229487	XS (46–50 cm)	●
NS229488	S (50–54 cm)	●
NS219410	M (54–58 cm)	●
NS229483	L (58–62 cm)	●

KIT CONTAINS:

- Electrode cap for 19-channel EEG recording (Electro-Cap)
- Body harness
- Electrode gel Unimax (250 g)
- Special needle for gel

Electrode cap Electro-Cap for 19-channel EEG recording



ORDER CODE	SIZE	COLOUR
NS103130	XS (46–50 cm)	●
NS102949	S (50–54 cm)	●
NS102383	M (54–58 cm)	●
NS102950	L (58–62 cm)	●

Ear electrode E5-9S for Electro-Cap



ORDER CODE: NS103620

Electrode system MCScap-EC21



ORDER CODE: NS226428

KIT CONTAINS:

- Electrode system MCScap-EC21
- Electrode gel Unimax (250 g)
- Syringe (20 ml)
- Special needle for gel

SIZE:

- Inf II (26–30 cm)
- Inf I (32–36 cm)
- XS (36–42 cm)
- S (42–48 cm)
- M (48–54 cm)
- L (54–60 cm)

Electrode system MCScap-NS5-2515



ORDER CODE: NS233522

KIT CONTAINS:

- Electrode system MCScap-NS5-2515
- Electrode gel Unimax (250 g)
- Syringe (20 ml)
- Special needle for gel

SIZE:

- S (42–48 cm)
- M/S (45–51 cm)
- M (48–54 cm)
- L/M (51–57 cm)
- L (54–60 cm)
- XL/L (57–63 cm)
- XL (60–66 cm)

Electrode cap MCScap-52



ORDER CODE: NS222647

SIZE:

- M (48–54 cm)
- L (54–60 cm)

Electrode system MCScap-52



ORDER CODE: NS104019

KIT CONTAINS:

- Textile cap with fixing rings (modified pattern 10-10) — 2 pcs.
- EEG electrode for MCScap, Ag/AgCl — 48 pcs.
- Ground EEG electrode for MCScap, Ag/AgCl
- Reserve EEG electrode for MCScap, Ag/AgCl — 3 pcs.
- Ear electrode MCScap-A — 2 pcs.
- Electrode gel Unimax (250 g)
- Syringe (20 ml)
- Special needle for gel

SIZE:

- S (42–48 cm) + M (48–54 cm)
- M (48–54 cm) + L (54–60 cm)

Adapter for connecting amplifiers to MCScap electrodes via integrated connector



ORDER CODE: NS220431

EEG electrode for MCScap



ORDER CODE: NS223256

Ear adapter MCScap



ORDER CODE: NS223357

Elastic bandage for electrode cap fixation



ORDER CODE: NS990998.027

QTY: set of 3 pcs.

VIDEO EEG ACCESSORIES

Built-in MIC IR 4x zoom dome camera on a tripod



ORDER CODE: NS223395

KIT CONTAINS:

- PTZ IP Dahua 4x zoom dome camera (NS233526)
- Video camera tripod TV-300A (NS226058)
- Splitter cable (NS015103.036)
- Mains cable
- Extension cable (NS015103.035)
- Clamp plate for camera tripod (NS007998.004)
- "PTZ IP Dahua dome camera" guidelines

PTZ IP Dahua 4x zoom dome camera



ORDER CODE: NS233526

Reversing coloured video camera with PSU



ORDER CODE: NS223396

KIT CONTAINS:


- Video camera Sony SNC-EP50 (NS228480)
- Power unit D-Link (NS228481)

Video capture device



ORDER CODE: NS201556

EEG cap extension cable



ORDER CODE: NS007103.023-30

LENGTH: 3 m

USB extension cable, A-B



ORDER CODE: NS215431

LENGTH: 20 m

Adapter USB→COM TU-S9



ORDER CODE: NS225523

LENGTH: 0.6 m

Power unit for microphone and IR projector



ORDER CODE: NS223557

STIMULATORS

LED photic stimulator PhS-1



ORDER CODE: NS005302.005

LED photic stimulator PhS-2




ORDER CODE: NS066201.002

Visual stimulator (LED goggles)



ORDER CODE: NS005302.001

Adapter for high resolution pattern-stimulator with dual monitor video card



ORDER CODE: NS223949

KIT CONTAINS:

- Dual monitor video card (NS208969)
- Adapter for high resolution pattern-stimulator (NS033201.003)

Adapter for pattern-stimulator connection



ORDER CODE: NS033201.005

SVGA extension cable



ORDER CODE: NS201808

LENGTH: 3 m

Auditory stimulator TDH-39 (headphones)



ORDER CODE: NS032305.010

CONNECTOR: 3.5 stereo

Auditory stimulator TA-01 (headphones)




ORDER CODE: NS015305.001

CONNECTOR: 3.5 stereo

MISCELLANEOUS

The coloured trim for Neuron-Spectrum-4/P front panel



ORDER CODE: NS015210.035-06

Patient button




ORDER CODE	CONNECTOR
NS028201.006	USB
NS034201.003	2.5 stereo

Trigger unit TU-1



ORDER CODE: NS990201.006

Fixing ring for bridge EEG electrode



ORDER CODE: NS224634

CFM pod



ORDER CODE: NS015201.040-02.20

LENGTH: 2 m

USB hub with cable KM-7-2



ORDER CODE: NS042999.002-02

KIT CONTAINS:

- USB hub KM-7-2 (NS042201.009)
- USB cable (1.8 m) (NS007103.005-01)
- Power supply cord (3 m)
- Technical manual

EMG ACCESSORIES



BASIC ACCESSORIES

Cable for disposable NCS electrodes



ORDER CODE: NS232891

KIT CONTAINS:

- Cable for surface, ring EMG electrodes (1.5 m) (NS990103.003-15)
- Cable for disposable electrode Alligator clip → touch-proof (0.2 m) — 2 pcs. (NS990103.027-03.02, NS990103.027-04.02)

Surface electrode with cable SE-1



ORDER CODE: NS990998.001

KIT CONTAINS:

- Cable for surface, ring EMG electrodes (1.5 m) (NS990103.003-15)
- Cup electrode for surface electrode SE-1 (0.2m) — 2 pcs. (NS990106.008-03.02, NS990106.008-04.02)

Surface electrode SE-1



ORDER CODE	COLOUR
NS990106.008-03.02	●
NS990106.008-04.02	●

LENGTH: 0.2 m

Ring electrode (narrow) with cable RE-1



ORDER CODE: NS990998.005

KIT CONTAINS:

- Cable for surface, ring EMG electrodes (0.5 m) (NS990103.003-05)
- Ring electrode (narrow) RE-1 (0.2 m) — 2 pcs. (NS990106.015-01, NS990106.015-02)

Ring electrode (narrow) RE-1



ORDER CODE	COLOUR
NS990106.015-01	●
NS990106.015-02	●

LENGTH: 0.2 m

Ring electrode (wide) with cable RE-2



ORDER CODE: NS990998.004

KIT CONTAINS:

- Cable for surface, ring EMG electrodes (1.5 m) (NS990103.003-15)
- Ring electrode (wide) RE-2 (0.2 m) — 2 pcs. (NS990106.007-01, NS990106.007-02)

Ring electrode (wide) RE-2



ORDER CODE	COLOUR
NS990106.007-01	●
NS990106.007-02	●

Cable for surface, ring EMG electrodes



ORDER CODE: NS990103.003-15

LENGTH: 1.5 m

Bar electrode (adult) BE-2



ORDER CODE: NS990998.029

KIT CONTAINS:

- Bar electrode (adult) BE-2
- Velcro strap (300 mm) (NS990211.002)

Bar electrode (pediatric) BE-1



ORDER CODE: NS990998.028

KIT CONTAINS:

- Bar electrode (pediatric) BE-1
- Velcro strap (300 mm) (NS990211.002)

Cable for disposable electrode ("Alligator" clip → touch-proof)



ORDER CODE	LENGTH	COLOUR
NS990103.027-03.02	0.2 m	●
NS990103.027-04.02	0.2 m	●
NS990103.027-01.10	1 m	○
NS990103.027-02.10	1 m	●
NS990103.027-03.10	1 m	●
NS990103.027-04.10	1 m	●
NS990103.027-05.10	1 m	●
NS990103.027-06.10	1 m	●
NS990103.027-02.12	1.2 m	●
NS990103.027-03.12	1.2 m	●
NS990103.027-04.12	1.2 m	●

Reusable concentric needle electrode D1216



ORDER CODE	NEEDLE	COLOUR
NS990106.039-001	20 mm/0.4 mm	●
NS990106.039-003	30 mm/0.4 mm	●
NS990106.039-004	40 mm/0.5 mm	●
NS990106.039-005	50 mm/0.45 mm	●
NS990106.039-006	60 mm/0.6 mm	●

Disposable concentric needle electrode B50600



ORDER CODE	NEEDLE	COLOUR
NS990106.038-001	25 mm/0.3 mm	●
NS990106.038-002	25 mm/0.4 mm	●
NS990106.038-003	37 mm/0.45 mm	●
NS990106.038-004	50 mm/0.45 mm	●
NS990106.038-005	75 mm/0.6 mm	●

Disposable hypodermic needle electrode with cable L2575



ORDER CODE	NEEDLE	COLOUR
NS990106.040-001	25 mm/0.3 mm	●
NS990106.040-002	25 mm/0.4 mm	●
NS990106.040-003	37 mm/0.4 mm	●
NS990106.040-004	37 mm/0.45 mm	●
NS990106.040-005	50 mm/0.5 mm	●
NS990106.040-006	75 mm/0.7 mm	●

Reusable cup electrode with cable C1234 (silver)



ORDER CODE	DIAMETER	LENGTH
NS990106.049-001	6 mm	1 m
NS990106.049-002	6 mm	1.5 m
NS990106.049-003	6 mm	2 m
NS990106.049-004	10 mm	1 m
NS990106.049-005	10 mm	1.5 m
NS990106.049-006	10 mm	2 m

- OTHER CUP COATING:**
- silver-chloride
 - gold

Cup electrode with cable



ORDER CODE: NS990106.028

LENGTH: 1 m

CAN BE SUPPLIED AS SET:

- 5 pcs., colours: ● ● ● ● ○ (NS990998.023)
- 8 pcs., colours: ● ● ● ● ○ ● ● ● (NS990998.024)

Reusable cable for concentric needle electrode D60451



ORDER CODE: NS990103.044

Disposable surface electrode F3001



ORDER CODE: NS104335

SIZE: 22x34 mm

QTY: 100 pcs.

Ground electrode with cable



ORDER CODE	MODEL	SIZE
NS990998.007	GE-1	250 mm
NS990998.006	GE-2	400 mm
NS990998.015	GE-3	700 mm

Y-adaptor



ORDER CODE: NS006103.019

Cable for OAE probe connection



ORDER CODE: NS006103.024

Adapter DIN → touch-proof



ORDER CODE: NS065221.005

STIMULATORS

Adapter for high resolution pattern-stimulator with dual monitor video card



ORDER CODE: NS223949

KIT CONTAINS:

- Dual monitor video card (NS208969)
- Adapter for high resolution pattern-stimulator (NS033201.003)

Adapter for pattern-stimulator connection



ORDER CODE: NS033201.005

Visual stimulator (LED goggles)



ORDER CODE: NS005302.001

SVGA extension cable



ORDER CODE: NS201808

LENGTH: 3 m

Auditory stimulator TDH-39
(headphones)



ORDER CODE: NS032305.010

CONNECTOR: 3.5 stereo

Auditory stimulator TDH-39
(headphones)



ORDER CODE: NS032305.009

CONNECTOR: 2x6.3 mono

Auditory stimulator TDH-39
(headphones, shielded)



ORDER CODE: NS032305.012

CONNECTOR: 3.5 stereo

Auditory stimulator TDH-39
(headphones, shielded)



ORDER CODE: NS032305.011

CONNECTOR: 2x6.3 mono

Auditory stimulator TA-01
(headphones)



ORDER CODE: NS015305.001

CONNECTOR: 3.5 stereo

Auditory stimulator HDA-280
(headphones)



ORDER CODE: NS032305.003

CONNECTOR: 3.5 stereo

Insert earphones ER3C-10



ORDER CODE: NS229864

CONNECTOR: 2x6.3 mono

KIT CONTAINS:

- Insert earphones ER3-10
- Cable for insert earphones
- Disposable foam eartip (10 mm) — 20 pcs.
- Disposable foam eartip (13 mm) — 20 pcs.
- Disposable foam eartip (18 mm) — 2 pcs.
- Velcro strap
- "Insert Earphones ER3-10" guidelines

Adapter for insert earphones
(3.5 stereo → 2x6.3 mono)



ORDER CODE: NS032103.004

Set of eartips



ORDER CODE	SIZE
NS007998.002	adult
NS007998.001	pediatric

OAE probe OAE-02 (mini-DIN)



ORDER CODE: NS006355.003-01

CONNECTOR: mini-DIN

OAE probe tip



ORDER CODE: NS006221.001

Test cavity



ORDER CODE: NS006201.013

Adjustable electro stimulating probe
AESP-1



ORDER CODE: NS056998.001

Adjustable electro stimulating probe
AESP-2



ORDER CODE: NS990106.037

Stimulation electrode with steel
stimulation point (pediatric) SSE-1



ORDER CODE: NS990106.076-025

Stimulation electrode with steel
stimulation point (adult) SSE-2



ORDER CODE: NS990106.075-025

Stimulating bar electrode (adult) SBE-2



ORDER CODE: NS990998.009

Stimulating bar electrode with
replaceable steel and felt stimulation
pads (adult) SBE-2



ORDER CODE: NS990998.035

Stimulating bar electrode with
replaceable steel and felt stimulation
pads (pediatric) SBE-1



ORDER CODE: NS990998.008

Tendon hammer



ORDER CODE: NS040356.001

CONNECTOR: USB

Temperature sensor



ORDER CODE: NS039351.003

CONNECTOR: USB

Neuro-MEP auditory-visual stimulator unit



ORDER CODE: NS025201.013

Neuro-MEP electrical stimulator control unit



ORDER CODE: NS024201.017

MISCELLANEOUS

Dedicated keyboard DK-01



ORDER CODE: NS103966

KIT CONTAINS:

- Dedicated keyboard DK-01 (NS035201.005)
- USB — Bluetooth adapter (NS217287)
- A-B USB cable (1.8 m) (NS007103.005-01)
- Technical manual "DK-01"

Clamp for dedicated keyboard



ORDER CODE: NS035221.001

Dedicated keyboard DK-02



ORDER CODE: NS232332

KIT CONTAINS:

- Dedicated keyboard DK-02 (NS035201.006)
- USB — Bluetooth adapter (NS217287)
- USB cable (1.8 m) (NS007103.005-01)
- Technical manual "DK-02"

Footswitch



ORDER CODE	CONNECTOR
NS028353.004	USB
NS034353.003	3.5 mm stereo

Patient button



ORDER CODE	CONNECTOR
NS028201.006	USB
NS034201.003	2.5 stereo

USB hub with cable KM-7-2



ORDER CODE: NS042999.002-02

KIT CONTAINS:


- USB hub KM-7-2 (NS042201.009)
- USB cable (1.8 m) (NS007103.005-01)
- Power supply cord (3 m)
- Technical manual

Trigger unit TU-1



ORDER CODE: NS990201.006


Reusable limb clamp ECG electrode



ORDER CODE	MODEL	SIZE
NS103212	F9024SSC	adult
NS101823	F9023	pediatric

QTY: 4 pcs.

3-lead ECG cable



ORDER CODE: NS007103.016

Measuring reel



ORDER CODE: NS224422

Marker pen



ORDER CODE: NS224039

Medical tape Transpore



ORDER CODE: NS228281

SIZE: 5 cm x 9.1 m

IONM ACCESSORIES




BASIC ACCESSORIES

Neuro-TES transcranial electrical stimulator



ORDER CODE: NS070201.004-020

Patient simulator



ORDER CODE: NS233345

KIT CONTAINS:

- Patient simulator (NS064201.011-02)
- Cable to attach transcranial electrical stimulator switch to patient simulator (NS064103.010)
- Cable to attach acquisition pod to patient simulator (NS064103.011)
- Cable to attach electrical stimulator pod to patient simulator (NS064103.012)
- Touch-proof — touch-proof cable (NS990103.038-03-05)

Neuro-IOM 9-channel electrical stimulator




ORDER CODE: NS061201.011-02

Neuro-MEP low current electrical stimulator



ORDER CODE: NS024201.019

Neuro-TES electronic switch



ORDER CODE: NS070201.007-02.50

Neuro-IOM low current stimulator pod



ORDER CODE: NS061201.009-02.50

Neuro-IOM electrical stimulator pod



ORDER CODE	TYPE	COLOUR
NS061201.007-12	for 1-4 channels	●
NS061201.007-22	for 5-8 channels	●

Neuro-IOM acquisition pod



ORDER CODE	MODEL	COLOUR
NS064201.009-02	A (1–8), 5 m	●
NS064201.009-12	B (9–16), 5 m	●
NS064201.009-22	C (17–24), 5 m	●
NS064201.009-32	D (25–32), 5 m	●

MMG unit



ORDER CODE: NS077201.001-02

MMG sensor



ORDER CODE: NS077356.001

Adapter



ORDER CODE	MODEL
NS990103.052	4-to-1 adapter
NS990103.057	8-to-1 adapter

ES detector



ORDER CODE: NS074999.001

- KIT CONTAINS:**
- ES detector (NS074201.001-02)
 - ES detector clip — 2 pcs. (NS074354.001-50)

ES detector clip



ORDER CODE: NS074354.001-50

ELECTRODES

Disposable subdermal corkscrew needle electrode S50715



ORDER CODE: NS990106.041-001

LENGTH: 1 m

Disposable subdermal curved needle electrode S45638



ORDER CODE: NS990106.045-001

LENGTH: 1 m

Disposable subdermal hooked needle electrode S50719



ORDER CODE: NS990106.047-001

LENGTH: 1 m

Disposable subdermal single needle electrode S50716



ORDER CODE: NS990106.043-001

LENGTH: 1 m


Disposable subdermal single needle electrode S50716 (twisted pair)



ORDER CODE: NS990106.043-06

LENGTH: 1 m

Disposable concentric probe D3600



ORDER CODE: NS990106.052

SIZE: 45 mm

Disposable monopolar probe D3602 (45 mm)



ORDER CODE: NS990106.054

SIZE: 45 mm


Disposable pedicle screw probe D3603



ORDER CODE: NS990106.055

SIZE: 45 mm

Disposable bipolar probe D3604



ORDER CODE: NS990106.056

SIZE: 45 mm


Disposable monopolar probe, tapered D3605



ORDER CODE: NS990106.057

SIZE: 45 mm

Disposable bipolar laryngeal electrode



ORDER CODE: NS990106.087

Disposable IOM subdural electrode



ORDER CODE	MODEL	LENGTH
NS220011	1-strip 4-contacts	0.5 m
NS227557	1-strip 8-contacts	1.2 m
NS227556	2-strip 16-contacts	1.2 m
NS228811	4-strip 16-contacts	1.2 m

Vagus stimulation electrode



ORDER CODE: NS064106.001

IOM probe handle with controls



ORDER CODE: NS061103.001

PROBES:

- Ball-tip stimulation probe (NS233517)
- Disposable monopolar stimulation probe (NS233516)

Ground electrode with cable



ORDER CODE	MODEL	SIZE
NS990998.007	GE-1	25 cm
NS990998.006	GE-2	40 cm
NS990998.015	GE-3	70 cm

Cable for stimulating electrodes with touch-proof connector



ORDER CODE: NS070103.002-045

LENGTH: 4.5 m

Insert earphones ER3C-10



ORDER CODE: NS229864

KIT CONTAINS:

- Insert earphones ER3-10
- Cable for insert earphones
- Disposable foam eartip (10 mm) — 20 pcs.
- Disposable foam eartip (13 mm) — 20 pcs.
- Disposable foam eartip (18 mm) — 2 pcs.
- Velcro strap
- “Insert Earphones ER3-10” guidelines

CONNECTOR: 2x6.3 mono

Adapter for insert earphones (3.5 stereo → 2x6.3 mono)



ORDER CODE: NS032103.004

Visual stimulator PhS-4 (LED goggles)



ORDER CODE: NS064302.001

MISCELLANEOUS

Stereo audio cable (2RCA → 2RCA)



ORDER CODE: NS213199

Video cable



ORDER CODE: NS015103.013

LENGTH: 13 m

Measuring reel



ORDER CODE: NS224422

Medical tape Transpore



ORDER CODE: NS228281

SIZE: 5 cm × 9.1 m

ABR AND OAE ACCESSORIES



BASIC ACCESSORIES

Adapter for 2-channel EP recording



ORDER CODE: NS032103.006-05

Neuro-Audio-Screen cable for EP electrodes



ORDER CODE: NS057103.008-05

Y-adapter



ORDER CODE: NS006103.019

ECochG electrode Tiptrode



ORDER CODE	MODEL	SIZE
NS104334	ER3-26B	10 mm
NS103769	ER3-26A	13 mm

Disposable surface electrode F3001



ORDER CODE: NS104335

SIZE: 22x34 mm
QTY: 100 pcs.

Cable for disposable electrode (Alligator clip → touch-proof)



ORDER CODE	LENGTH	COLOUR
NS990103.027-03.02	0.2 m	●
NS990103.027-04.02	0.2 m	●
NS990103.027-02.05	0.5 m	●
NS990103.027-03.05	0.5 m	●
NS990103.027-04.05	0.5 m	●
NS990103.027-06.05	0.5 m	●
NS990103.027-01.10	1 m	○
NS990103.027-02.10	1 m	●
NS990103.027-03.10	1 m	●
NS990103.027-04.10	1 m	●
NS990103.027-05.10	1 m	●
NS990103.027-06.10	1 m	●
NS990103.027-02.12	1.2 m	●
NS990103.027-03.12	1.2 m	●
NS990103.027-04.12	1.2 m	●

Cable for disposable electrode connection (button clip → touch-proof)



ORDER CODE	LENGTH	COLOUR
NS990103.026-01.05	0.5 m	○
NS990103.026-02.05	0.5 m	●
NS990103.026-03.05	0.5 m	●
NS990103.026-04.05	0.5 m	●
NS990103.026-06.05	0.5 m	●
NS990103.026-01.10	1 m	○
NS990103.026-02.10	1 m	●
NS990103.026-03.10	1 m	●
NS990103.026-04.10	1 m	●
NS990103.026-06.10	1 m	●

Cable for ECochG electrode



ORDER CODE: NS032103.005-04.10

Cup electrode with cable



ORDER CODE: NS990106.028

LENGTH: 1 m

CAN BE SUPPLIED AS SET:

- 5 pcs., colours: ● ● ● ● ○ (NS990998.023)
- 8 pcs., colours: ● ● ● ● ● ● ● ● (NS990998.024)

Reusable cup electrode with cable C1234 (silver)



ORDER CODE	DIAMETER	LENGTH
NS990106.049-001	6 mm	1 m
NS990106.049-002	6 mm	1.5 m
NS990106.049-003	6 mm	2 m
NS990106.049-004	10 mm	1 m
NS990106.049-005	10 mm	1.5 m
NS990106.049-006	10 mm	2 m

OTHER CUP COATING:

- silver-chloride
- gold

STIMULATORS

Auditory stimulator TDH-39 (headphones)



ORDER CODE: NS032305.009

CONNECTOR: 2x6.3 mono

Auditory stimulator HDA-280 (headphones)



ORDER CODE: NS032305.003

CONNECTOR: 3.5 stereo

Auditory stimulator TDH-39 (headphones, shielded)



ORDER CODE: NS032305.011

CONNECTOR: 2x6.3 mono

Auditory stimulator TA-01 (headphones)



ORDER CODE: NS015305.001

CONNECTOR: 3.5 stereo

Bone vibrator



ORDER CODE: NS032998.001

CONNECTOR: 6.3 mono

Insert earphones ER3C-10



ORDER CODE: NS229864

CONNECTOR: 2x6.3 mono

KIT CONTAINS:

- Insert earphones ER3-10
- Cable for insert earphones
- Disposable foam eartip (10 mm) — 20 pcs.
- Disposable foam eartip (13 mm) — 20 pcs.
- Disposable foam eartip (18 mm) — 2 pcs.
- Velcro strap
- "Insert Earphones ER3-10" guidelines

Adapter for insert earphones (3.5 stereo → 2x6.3 mono)



ORDER CODE: NS032103.004

OAE probe OAE-02



ORDER CODE: NS006355.003-01

CONNECTOR: mini-DIN

OAE probe OAE-05



ORDER CODE: NS075355.003-04

CONNECTOR: micro-HDMI

OAE probe OAE-04-2 (touchscreen)



ORDER CODE: NS057355.002-00

CONNECTOR: binder

OAE probe OAE-03-2 (touchscreen)



ORDER CODE: NS057355.001-00

CONNECTOR: mini-DIN

OAE probe tip



ORDER CODE: NS006221.001

Set of eartips



ORDER CODE	MODEL
NS007998.002	universal
NS007998.001	pediatric

Eartip ER3-14B (pediatric)



ORDER CODE: NS104278

SIZE: 10 mm

Eartip ER3-14A (adult)



ORDER CODE: NS104277

SIZE: 13 mm

Eartip



ORDER CODE	MODEL	SIZE
NS222834	KR004	4 mm
NS222835	KR005	5 mm
NS221417	KR006	6 mm
NS222836	KR007	7 mm
NS221418	KR008	8 mm
NS221419	KR009	9 mm
NS222837	KR010	10 mm
NS221420	KR011	11 mm
NS221421	KR012	12 mm
NS221422	KR013	13 mm
NS221423	KR014	14 mm
NS222875	KR016	16 mm
NS221416	KR036	3–6 mm
NS222833	MF003	3–8 mm

MISCELLANEOUS

Test cavity



ORDER CODE: NS006201.013

Patient button



ORDER CODE	CONNECTOR
NS028201.006	USB
NS034201.003	2.5 stereo

Thermal printer with wireless Bluetooth interface, Custom



ORDER CODE: NS228532

Thermal printer with wireless Bluetooth interface, Citizen



ORDER CODE: NS222996

Power supply unit



ORDER CODE: NS057201.009

Loopback cable



ORDER CODE: NS032103.012

USB — Bluetooth adapter



ORDER CODE: NS217287

Probe tip extractor



ORDER CODE: NS006206.016

Dental floss Oral-B Superfloss



ORDER CODE: NS215488

QTY: 50 pcs.

ERG ACCESSORIES



BASIC ACCESSORIES

Corneal loop electrode



ORDER CODE: NS006106.005

Corneal hook electrode



ORDER CODE: NS006106.004

Adapter for corneal electrode



ORDER CODE: NS006103.011

Ear EEG electrode



ORDER CODE: NS015106.015

Cup electrode with cable



ORDER CODE: NS990106.028

LENGTH: 1 m

CAN BE SUPPLIED AS SET:

- 5 pcs., colours: ● ● ● ● ○ (NS990998.023)
- 8 pcs., colours: ● ● ● ● ● ● ● ● (NS990998.024)

Y-adapter



ORDER CODE: NS006103.019

Reusable cup electrode with cable C1234 (silver)



ORDER CODE	DIAMETER	LENGTH
NS990106.049-001	6 mm	1 m
NS990106.049-002	6 mm	1.5 m
NS990106.049-003	6 mm	2 m
NS990106.049-004	10 mm	1 m
NS990106.049-005	10 mm	1.5 m
NS990106.049-006	10 mm	2 m

OTHER CUP COATING:

- silver-chloride
- gold

Dedicated keyboard DK-02



ORDER CODE: NS232332

KIT CONTAINS:

- Dedicated keyboard DK-02 (NS035201.006)
- USB — Bluetooth adapter (NS217287)
- USB cable (1.8 m) (NS007103.005-01)
- Technical manual "DK-02"

Holder for penlights and mini-ganzfeld stimulator



ORDER CODE: NS025201.012


STIMULATORS

Visual stimulator (LED goggles)



ORDER CODE: NS005302.001

Set of penlights with concentrators



ORDER CODE: NS006302.004

COLOUR:




Mini-ganzfeld stimulator



ORDER CODE: NS025302.001

Adapter for high resolution pattern-stimulator with dual monitor video card



ORDER CODE: NS223949

KIT CONTAINS:

- Dual monitor video card (NS208969)
- Adapter for high resolution pattern-stimulator (NS033201.003)

MISCELLANEOUS

Dedicated keyboard DK-01



ORDER CODE: NS103966

KIT CONTAINS:

- Dedicated keyboard DK-01 (NS035201.005)
- USB — Bluetooth adapter (NS217287)
- A-B USB cable (1.8 m) (NS007103.005-01)
- Technical manual "DK-01"

USB hub with cable KM-7-2



ORDER CODE: NS042999.002-02

KIT CONTAINS:

- USB hub KM-7-2 (NS042201.009)
- USB cable (1.8 m) (NS007103.005-01)
- Power supply cord (3 m)
- Technical manual

SVGA extension cable



ORDER CODE: NS201808

LENGTH: 3 m

Holder for penlights and mini-ganzfeld stimulator



ORDER CODE: NS025201.012

ACCESSORIES
FOR MAGNETIC
STIMULATORS



COILS FOR MAGNETIC STIMULATORS



Small ring
100 mm



Big ring
150 mm



Small figure-of-eight
50 mm



Figure-of-eight
100 mm



Angulated figure-of-eight
100 mm



Double cone
125 mm



Non-cooled

RC-02-100
(NS058304.001)

RC-02-150
(NS058304.002)

FEC-02-50
(NS058304.005)

FEC-02-100
(NS058304.013)

FEC-02-100-P placebo
(NS058304.004)

AFEC-02-100
(NS058304.014)

AFEC-02-100-P placebo
(NS058304.009)

DCC-02-125
(NS058304.017)



Cooled

RC-02-150-C
(NS058304.003)

FEC-02-100-C
(NS058304.010)

AFEC-02-100-C
(NS058304.012)

AFEC-02-100-C-P placebo
(NS058304.011)

DCC-03-125-C
(NS063304.016)

BASIC ACCESSORIES

Control cable for Neuro-MS/D cooling unit (DB-15M → DB-9F)



ORDER CODE: NS058103.070

Control cable for Neuro-MS/D cooling and expansion units



ORDER CODE: NS058103.071

HV cable for Neuro-MS/D expansion unit connection



ORDER CODE: NS058103.101

HV cable for Neuro-MS/D extra power supply unit connection



ORDER CODE: NS058103.102

High-voltage cable for Neuro-MS extra unit connection



ORDER CODE: NS063103.007

End cap for Neuro-MS/D



ORDER CODE: NS058105.002

End cap for Neuro-MS




ORDER CODE: NS063105.001

Trig-in cable for synchronization via electrical stimulator (DIN5 240° → DIN8)



ORDER CODE: NS014103.032

Trigger unit TU-2



ORDER CODE: NS990201.007-02

Equipotential cable



ORDER CODE: NS058103.032-050

LENGTH: 0.5 m

Hydraulic system adapter



ORDER CODE: NS058214.001

Hydraulic connecting hose



ORDER CODE: NS230521

Connecting pipe



ORDER CODE: NS230207

Trig-in cable



ORDER CODE	SYNCHRONIZATION TYPE
NS014103.039	DIN8 ← BNC
NS014103.047	DIN8 → BNC

Footswitch for Neuro-MS



ORDER CODE: NS063353.001

Coil holder (trolley/wall mounted) K-5



ORDER CODE: NS058998.003

Coil holder (trolley/wall mounted) K-8



ORDER CODE: NS059998.004

Silicone oil



ORDER CODE: NS058998.004

VOLUME: 5 l

Patient cap



ORDER CODE	SIZE
NS058213.001	42-54
NS058213.001-01	54-66

ECG ACCESSORIES



BASIC ACCESSORIES

Disposable ECG electrode F9079



ORDER CODE: NS104391

SIZE: 43x45 mm
QTY: 50 pcs.

Disposable ECG electrode 10C



ORDER CODE: NS226828

DIAMETER: 26 mm
QTY: 30 pcs.

Adapter for disposable ECG electrode with button clip



ORDER CODE: NS102063

COLOUR:
● ● ● ○

Reusable limb clamp ECG electrode



ORDER CODE	MODEL	SIZE
NS103212	F9024SSC	adult
NS101823	F9023	pediatric

QTY: 4 pcs.

Reusable precordial suction chest ECG electrode



ORDER CODE	MODEL	SIZE
NS103200	F9009SSC	adult
NS101822	F9008	pediatric

DIAMETER: 24 mm
MATERIAL: Ag/AgCl

Reusable peripheral plate electrode for limb



ORDER CODE: NS103642

SIZE: 30x20 mm

Reusable peripheral plate electrode for chest



ORDER CODE: NS102080

SIZE: 15x20 mm

Rubber band for ECG electrode fixation on limb



ORDER CODE: NS011211.002

SIZE: 300x25 mm

Rubber belt for ECG electrode fixation on the chest



ORDER CODE: NS011211.007

Plastic button to fasten rubber belt PG-905/99



ORDER CODE: NS226453

Set of electrodes for exercise testing



ORDER CODE: NS101922

KIT CONTAINS:

- Reusable peripheral plate electrode for chest (15×20 mm) — 6 pcs. (NS102080)
- Reusable peripheral plate electrode for limb (30×20 mm) — 4 pcs. (NS103642)
- Belt for ECG electrode fixation on chest (six sets of holes, 100×1500 mm) (NS226452)
- Rubber tape for ECG electrode fixation on limb (300×25 mm) (NS011211.002)
- Plastic button to fasten rubber belt — 3 pcs. (NS226453)

Set of electrodes for ECG recording according to Frank lead system



ORDER CODE: NS102970

KIT CONTAINS:

- Reusable peripheral plate electrode for limb (30×20 mm) — 6 pcs. (NS103642)
- Belt for ECG electrode fixation on chest (six sets of holes, 100×1500 mm) (NS226452)
- Plastic button to fasten rubber belt — 3 pcs. (NS226453)

3-lead ECG cable



ORDER CODE: NS007103.016

Poly-Spectrum-12/E cable for additional ECG leads



ORDER CODE: NS222044

4-lead ECG cable with banana connector



ORDER CODE	TYPE	LENGTH
NS222048	for Poly-Spectrum-8/V	2 m
NS222046	for VNS-Rhythm, VNS-Micro, Poly-Spectrum-8/V	1.4 m

10-lead ECG cable with banana connector



ORDER CODE	TYPE	LENGTH
NS104056	for Poly-Spectrum-12/E	3.1 m
NS103272	for Poly-Spectrum-8, 8/E, 8/EX	1.3 m

Airflow sensor AS-4-20 for Poly-Spectrum-8, 8/E, VNS-Rhythm, VNS-Micro



ORDER CODE: NS990351.005

CONNECTOR: 3.5 stereo
LENGTH: 2 m

Airflow sensor for Poly-Spectrum-12/E, 8/EX



ORDER CODE	MODEL	LENGTH
NS990351.004	BS-3-20	2 m
NS990351.004-01	BS-3-04	0.4 m
NS990351.004-02	BS-3-12	1.2 m

CONNECTOR: 2.5 stereo

MISCELLANEOUS

Tapuz ECG Belt



ORDER CODE: NS990106.082

Cable for Kettler cycle ergometer connection to PC



ORDER CODE: NS007103.001

Cable for Ergoline cycle ergometer connection to PC



ORDER CODE: NS007103.002

Battery charger GP



ORDER CODE: NS217175

ACCESSORIES FOR SPIROMETERS



BASIC ACCESSORIES

Assembled flow converter



ORDER CODE: NS019201.011

Flange



ORDER CODE: NS019205.003

Semiring clamp



ORDER CODE: NS019205.004

Replaceable grid




ORDER CODE: NS019205.005

Rubber ring gasket



ORDER CODE: NS225215

Pulse tube



ORDER CODE: NS019201.006

Mouthpiece



ORDER CODE: NS019204.002

Nose clip



ORDER CODE: NS220183

Disposable cardboard mouthpiece



ORDER CODE: NS102851

Calibration pump



ORDER CODE: NS019201.007

VOLUME: 3 l

STANDS, ARMS, HOLDERS
AND TROLLEYS
PASTES AND GELS



STANDS, ARMS, HOLDERS

Holder H-1S



ORDER CODE: NS016201.038

Holder for LED photic stimulator



ORDER CODE: NS016221.016-01

Holder for penlights and mini-ganzfeld stimulator



ORDER CODE: NS025201.012

Small holder SH-1S



ORDER CODE: NS016201.039

Holder for Skybox



ORDER CODE: NS065221.004

Stand for electronic unit SU-3



ORDER CODE: NS016998.007

Stand for electronic unit SU-7



ORDER CODE: NS016998.013

Stand for LED photic stimulator SU-8



ORDER CODE: NS016201.042

Stand for electronic unit (for Neuro-MEP-5) SU-12



ORDER CODE: NS016998.016

Stand for top fixation of ECG cable



ORDER CODE: NS016998.005

Coil holder (trolley/wall mounted) K-5



ORDER CODE: NS058998.003

Coil holder (trolley/wall mounted) K-8



ORDER CODE: NS059998.004

Slider for Neuro-MEP-Micro



ORDER CODE: NS034221.001

Mounting bracket KK-1



ORDER CODE: NS016221.032

Electrode holder



ORDER CODE: NS016221.023

Flexible arm for coil positioning K-3



ORDER CODE: NS016998.012

Video camera tripod TV-300A



ORDER CODE: NS226058

Forehead and chin support assembly



ORDER CODE: NS016998.004

Mounting plate for two amplifier units



ORDER CODE: NS006200.003

Holder for Neuro-Audio-Screen



ORDER CODE: NS075221.001

Holder for aScreen



ORDER CODE: NS075211.001

Mount to fix the electronic unit on stand



ORDER CODE: NS016221.009

Hook to fix the electronic unit on holder



ORDER CODE: NS006200.002

TROLLEYS

Trolley T-3



ORDER CODE: NS059998.003

Trolley T-3 with holder



ORDER CODE: NS224245

- KIT CONTAINS:**
- Trolley T-3 (NS059998.003)
 - Holder H-4S (NS016201.044)

Trolley T-3 with holder and mounting bracket



ORDER CODE: NS224244

- KIT CONTAINS:**
- Trolley T-3 (NS059998.003)
 - Holder H-4S (NS016201.044)
 - Mounting bracket KK-1 — 4 pcs. (NS016221.032)

Trolley with custors for magnetic stimulator T-4



ORDER CODE: NS059998.006

Trolley T-5



ORDER CODE: NS059221.018

Trolley HighGradeTech



ORDER CODE: NS227146

Trolley with monitor mount HumanScale



ORDER CODE: NS228842

Chair Greiner



ORDER CODE: NS228837

PASTES AND GELS

Electrode gel Unimax



ORDER CODE: NS223722

WEIGHT: 250 g

Conductive electrode spray Unisprey



ORDER CODE: NS225064

WEIGHT: 200 g

Abrasive paste for skin preparation
EVERI



ORDER CODE: NS103717

WEIGHT: 160 g

Electrode adhesive paste TEN20



ORDER CODE: NS228002

WEIGHT: 114 g

CONTACTS



Neurosoft
5, Voronin str., Ivanovo, 153032, Russia
neurosoft.com



Sales
Phone: +7 4932 24-04-34
Fax: +7 4932 24-04-35
E-mail: info@neurosoft.com



Technical Support
Phone: +7 4932 59-21-12
Fax: +7 4932 24-04-37
E-mail: help@neurosoft.com
Skype: [neurosoft.service.centre](https://www.skype.com/join/neurosoft.service.centre)